

List of EUS-FNA needles successfully used* with AQ-Flex™ 19 N Confocal Miniprobes™

Pre-loading of the AQ-Flex™ 19 N inside the needle before the procedure is mandatory. Always use a locking device to lock the position of probe within needle.

PRODUCT NAME	BRAND	REF.	COMPAT	OPERATING PROCEDURE
EchoTip® Ultra Endoscopic Ultrasound Needle	Cook Medical	ECHO-19	✓	Most commonly used
Expect™ Needle 19 Flex & Expect™ Slimline (SL) 19G	Boston Scientific	5004 5553	✓	High flexibility helps for transduodenal approaches
SonoTip® Pro Control	MediGlobe GmbH	GUS-33-21-019	✓	Less commonly used
Beacon™ fine needle aspiration 19G	Medtronic	DSN-19-01	✓	Pre-load on white handle (needle)



* Only these needles have been tested for compatibility. If you want to use a needle that is not on the list, please ask your Mauna Kea Technologies contact or send an email to support@maunakeatech.com

List of EUS-FNA needles unsuccessfully used* with AQ-Flex™ 19 N Confocal Miniprobes™

PRODUCT NAME	BRAND	REF.	COMPAT	OPERATING PROCEDURE
EchoTip® Ultra Endoscopic Ultrasound Access Needle	Cook Medical	ECHO-HD-19-A	✗	 Not Recommended Pre-loading impossible (beveled stylet, flat needle)
EchoTip® ProCore	Cook Medical	ECHO-HD-19-C	✗	 Not Recommended Risk of damage with the distal hole in the needle
EZ Shot™ 3	Olympus	NA-U201H-8019	✗	 Not Recommended, Pre-loading may be hard or impossible.
Acquire® Flexible Needle	Boston Scientific	5558	✗	 Not Recommended Risk of damage with the bevel in the needle
19G FlexNeedle™ Clearview	CONMED	CV2119R	✗	 Not Recommended Pre-loading impossible

* Only these needles have been tested for compatibility. If you want to use a needle that is not on the list, please ask your Mauna Kea Technologies contact or send an email to support@maunakeatech.com

List of Aspirations needles successfully used* with AQ-Flex™ 19 N Confocal Miniprobe™ in pulmonary tract

Pre-loading of the AQ-Flex™ 19 N inside the needle before the procedure is mandatory. Always use a locking device to lock the position of probe within needle.

PRODUCT NAME	BRAND	REF.	COMPAT	OPERATING PROCEDURE
FlexNeedle® 18g	Broncus	10005-2	✔	

*** Only these needles have been tested for compatibility. If you want to use a needle that is not on the list, please ask your Mauna Kea Technologies contact or send an email to support@maunakeatech.com**

Cellvizio® I.V.E. Systems with Confocal Miniprobe™ are regulated Medical Device, CE marked (CE 0459) (Class IIa - NB : G-MED) and FDA cleared. Mention of products marketed by third party companies is for informational purpose only. Mauna Kea Technologies assumes no liability for the use and operation of these products. This list is provided for informational purpose only and is based on bench tests performed with AQ-Flex™ 19 Confocal Miniprobe™. It does not constitute an endorsement nor a recommendation and Mauna Kea Technologies does not guaranty the compatibility of the AQ-Flex™ 19 with these needles. Confocal Miniprobe™ availability cannot be guaranteed in all countries. For further information, please contact your local sales representative. The use of this medical device is exclusively reserved for healthcare professionals.

The Cellvizio® I.V.E. is a confocal laser system with fiber optic probes that are intended to allow imaging of the internal microstructure of tissues including, but not limited to, the identification of cells and vessels and their organization or architecture.

US Indications for use

Once connected to the Cellvizio® I.V.E. system, the AQ-Flex™ 19 N Confocal Miniprobe™ are intended to allow imaging of anatomical tracts, i.e., gastrointestinal and respiratory tracts, accessed by an endoscope, or endoscopic accessories (e.g. aspiration needles used during procedures including EUS-FNA, EBUS-TBNA and TBNA needles).

CE marked Indications for use

Once connected to the Cellvizio® I.V.E. system:

The AQ-Flex™ 19 N Confocal Miniprobe™ are intended to allow imaging of anatomical tracts, i.e., gastrointestinal and respiratory tracts, accessed by an endoscope, or endoscopic accessories, including through endoscopic needles.